



# **HiSPARC Station Problem and Solutions Documentation**

*Release 0.1 alpha*

**HiSPARC team**

November 29, 2012



# CONTENTS

|          |  |           |
|----------|--|-----------|
| <b>1</b> | <b>Introduction</b>                      | <b>3</b>  |
| <b>2</b> | <b>Frequently asked questions</b>        | <b>5</b>  |
| <b>3</b> | <b>Known issues</b>                      | <b>7</b>  |
| 3.1      | Structure . . . . .                      | 7         |
| <b>4</b> | <b>Software</b>                          | <b>9</b>  |
| 4.1      | HiSPARC Monitor does not start . . . . . | 9         |
| 4.2      | Hard Disc Space . . . . .                | 9         |
| 4.3      | HiSPARC DAQ Errors . . . . .             | 10        |
| 4.4      | Error in HiSPARC Monitor . . . . .       | 11        |
| <b>5</b> | <b>Hardware</b>                          | <b>13</b> |



The HiSPARC software can be downloaded here: [HiSPARC Installer](#)

A pdf version of this manual is available here: [HiSPARC maintenance manual](#)

For more information about HiSPARC, see the [HiSPARC Website](#).

For questions or suggestions please contact our project coordinator Surya Bonam ([info@hisparc.nl](mailto:info@hisparc.nl)).



# INTRODUCTION

This documentation describes the steps required to solve problems with an HiSPARC station. It first explains how to recognize different problems and then how to solve them.

The first part, *Frequently asked questions*, might help answer some questions that you have.

The second part, *Known issues*, describes some of the known issues that we have encountered with the HiSPARC software and hardware. It describes how to determine if the problem described is what you are experiencing and then possible solutions.





# FREQUENTLY ASKED QUESTIONS

- **Nothing works! What can I do?**
  - Don't Panic! Try to determine where the problem lies and find the solution in this documentation.
- **Where do I find the list of known possible issues?**
  - Here: *Known issues*
- **The histograms for my station don't look correct, how can I adjust the station settings to fix this?**
  - Read here about *calibrating the detectors*
- **I tried that, I didn't find a similar problem! What do I do now?**
  - Please contact your *cluster coordinator* for further assistance. They may be aware of such problems the corresponding solutions to the problem, and if they are not, they will consult us.



# KNOWN ISSUES

This is a list of known possible issues with HiSPARC stations. Each problem notes whether it has been fixed in a specific version of the software and if it was perhaps introduced by one. Also other parameters that could affect the occurrence of the issue are noted (e.g. Windows XP/7, RAM, cable lengths). Moreover, the Nagios warnings that an issue can trigger are mentioned.

---

**Note:** Multiple issues can cause the same Nagios warning.

---

## 3.1 Structure

Each problem has the following fields:

**First Sign** Explaining how you will probably notice the problem.

**Nagios** Nagios warnings that can be triggered.

**Determination** This is a small guide explaining how to make sure that the problem being described is what you are experiencing.

**Solution** How to solve it.

**Effects** The effects of this problem.



# SOFTWARE

This sections concerns itself with issues related to the HiSPARC station-software.

## 4.1 HiSPARC Monitor does not start

### 4.1.1 Missing directory

**First Sign** When the `STARThiSPARCSoftware` program runs and the other programs (HiSPARC DAQ and Updater) start normally but the HiSPARC Monitor does not appear or closes instantly.

**Nagios** `EventRate`, `StorageGrowth`, `StorageSize`, `TriggerRate`, and possibly `Buffer size Determination`

- Look in `hisparc/persistent/logs/src/` for the latest log file.
- Check if there is a line that contains the text `Error: unable to open database file.`
- Look in the `hisparc/persistent/data/` directory for a `hsmonitor` folder.
- If it does not exist than go to the solution, otherwise this is not the problem.

**Solution** Create the missing `hsmonitor` directory in `hisparc/persistent/data`.

**Effects** The missing directory causes the HiSPARC Monitor to be unable to store events in its SQLite database, preventing it from sending events to the Nikhef datastore. Note that the HiSPARC DAQ should be unaffected. No events should be lost, the DAQ will store events in its MySQL database until the hard disc fills up.

## 4.2 Hard Disc Space

### 4.2.1 To many logs

**First Sign** Nagios warning about Disc Space.

**Nagios** `Disc Space`.

**Determination**

- Look in hisparc/persistent/logs/.
- Check the size of the src directory by right-clicking on it and choosing ‘Properties’.
- Check if this is a significant fraction of the total disc space.

**Solution** Remove all logs from the src directory except for the current one (present date in dd-mm-yyyy.log). Select all (ctrl + a) logs in hisparc/persistent/logs/src. Deselect the current one (ctrl + click). Remove them using shift + delete (to bypass the Recycle Bin)

**Effects** If the disc is full the HiSPARC daq can not store events in the database, preventing the station from storing more events.

### 4.2.2 To many updaters

**First Sign** Nagios warning about Disc Space.

**Nagios Disc Space**

**Determination**

- Look in hisparc/persistent/downloads/.
- There should be some adminUpdater\_v##.zip and userUnpacker\_v##.exe files there.
- By right-clicking them you can see their file size is of the order of 100 MB.
- If there are many they can take up some space.

**Solution** Remove all userUnpacker and adminUpdater files except the newest ones. Do this by selecting them and pressing shift + delete to remove them directly.

**Effects** If the disc is full the HiSPARC daq can not store events in the database, preventing the station from storing more events.

## 4.3 HiSPARC DAQ Errors

### 4.3.1 Can not connect to buffer

**First Sign** Red LED in HiSPARC DAQ

**Nagios**

**Determination** From the Start menu start odbcad32.exe. Check if the hisparc buffer is there.

**Solution**

**Effects**

### 4.3.2 Not in DAQ Mode

**First Sign**

**Nagios TriggerRate**

**Determination** Look at the program HiSPARC DAQ, see if the button in the middle shows 'DAQ Mode'.

**Solution** Click the 'DAQ Mode' button in the HiSPARC DAQ.

**Effects** When the HiSPARC DAQ is not in DAQ Mode it will not store triggered events.

## 4.4 Error in HiSPARC Monitor

### 4.4.1 Malformed HisparcII.ini

**First Sign** Errors in the HiSPARC Monitor: `Uncatched exception in job:  
need more than 1 value to unpack. Restarting...`

**Nagios TriggerRate**

**Determination** Check in the file `hisparc/persistent/configuration/HisparcII.ini` for blank lines

**Solution** Remove any blank lines from `HisparcII.ini`

**Effects** Errors in the HiSPARC Monitor and no TriggerRate updates for Nagios.





# HARDWARE